


0.5 mL U-bottom Deep Well Storage Plate

MLS-DRWPU96-5



Description 0.5 mL U-bottom Deep Well Storage Plate, Non-Sterile

Materials PP (Polypropylene)

Features

- Autoclavable: no
- Can be stored at -196°C to +120°C
- Centrifuge to 4800 x G
- Alphanumeric grid-referencing rimmed version to stop cross contamination and enable a better seal
- Cylindrical well with conical bottom for optimal mixing and recovery
- No inner edges to allow better collection of magnetic beads
- Certified free of detectable RNase, DNase, DNA, PCR inhibitors, and tested pyrogen-free

| Item No. | Item | Sterile | Packaging Configuration |
|---------------|--|---------|-------------------------|
| MLS-DRWPU96-5 | 0.5 mL U-bottom Deep Well Storage Plate, 96 Round Well | No | 80/case |