

# 0.2 mL PCR Tube

## MLS-PCRT-02F



**Description** 0.2 mL Individual PCR Tube with Optically-Clear Flat Cap, Clear, Non-Sterile

**Materials** PP (Polypropylene)

- Features**
- Individual PCR tubes feature evaporation-resistant, easy-to-open-and-close lids that prevent sample loss
  - Autoclavable: yes
  - Ultra-thin, uniform wells ensure optimal heat transfer and reaction efficiency
  - Certified free of detectable RNase, DNase, DNA, PCR inhibitors, and tested pyrogen-free

Item No.	Item	Sterile	Packaging Configuration
MLS-PCRT-02F	0.2 mL Individual PCR Tube, Flat Cap	No	1000/unit 10 units/case

